



Dentium

White Seal



NEW

## Final prosthesis




# Final prosthesis

Measure depth



## White Seal

- Stably maintains shape, prevents top coating layer from collapsing
- Highly Biocompatible, made out of chemically stable Teflon material
- Less aesthetic concerns due to color change and odor
- One Length choice : 30 mm







# White Seal Application

Measure depth



Place the White Seal




# White Seal Application

Measure depth


Place the White Seal (2~3 mm downward form the crown top)



# White Seal Application



Easily cut



Useful for replenishment and removal



## White Seal Application [Clinical case]





## White Seal Application [Clinical case]

